


Invacare® Solace® Prevention 3080 Mattress


The soft heel section is upgraded to visco-elastic foam to cradle the sensitive area. This product provides comfort and support with built-in side supports and dual layer "fortress" cut foam.

Features

- Built-in side supports
- 70 Denier nylon top cover
- 11 oz. tri-laminate vinyl bottom
- Soft head and heel section
- Fortress cut convoluted foam



1. Test shows subject in supine position demonstrating contact measurement between key pressure points using FSA pressure mapping technology. Subject: Female, 5'2\", 105 lb.



2. Test shows subject in supine position demonstrating contact measurement between key pressure points using FSA pressure mapping technology. Subject: Female, 5'2\", 105 lb.



Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
SPS3080	80\" L 36\" W 6\" H	350 lb.	16 CFR 1632 16 CFR 1633	5 Years (Mattress) 1 Year (Cover)