

# Maple Leaf WHEELCHAIR



*The only limitation is your imagination*

**MLTR**<sup>500</sup>  
&  
**MLTR**<sup>600</sup>



TILT AND RECLINE  
MANUAL STEEL  
GERIATRIC WHEELCHAIRS



[www.mapleleafwheelchair.ca](http://www.mapleleafwheelchair.ca)

**MLTR<sup>500</sup>  
&  
MLTR<sup>600</sup>**

*The MLTR 500 is a geriatric tilt and recline steel frame manual wheelchair, lined with steel back panels for extra strength and comes with 4 casters, 3 total lock and 1 directional.*  
*The MLTR 600 differs in the fact that it includes 2 rear wheels instead of 2 casters in the back, with the frame designed to accommodate.*

## STANDARD FEATURES

### FRAME & UPHOLSTERY

- ♦ 16" and 18" chair widths
- ♦ 15" to 19" chair depths (infinite)
- ♦ MLTR 500
  - ♦ 19" and 21" seat heights (top of seat cushion)
- ♦ MLTR 600
  - ♦ 17.5", 18.5", 19.5" and 21" seat heights (top of seat cushion)
- ♦ 29" standard back height
- ♦ Durable chrome finish
- ♦ "Formfit" tray
- ♦ Unique one side adjustment for tray release, seat depth, legrest & footrest
- ♦ Adjustable, removable contoured headwings & headrest cushion
- ♦ Contoured upholstered seat & back with flip-down feature
- ♦ Vinyl upholstery in six colour choices

### TILT / RECLINE

- ♦ 0 - 30 degrees of infinite tilt adjustment
- ♦ 90 to 120 degrees of infinite recline adjustment
- ♦ Anti-tippers

### ARMRESTS & FOOTRESTS

- ♦ Durable padded arm pads
- ♦ Flip back removable arms
- ♦ Elevating legrest with fold-in footrest

### WHEELS

- ♦ MLTR 500
  - ♦ (3) - 5" total lock casters
  - ♦ (1) - 5" directional lock caster
- ♦ MLTR 600
  - ♦ 5", 6" or 8" casters
  - ♦ 20", 22" or 24" threaded urethane tires
- ♦ Durable chrome handrims



♦ ALL STEEL CONSTRUCTION

## ACCESSORIES & OPTIONS

- ♦ Non-standard seat width
- ♦ Non-standard back height
- ♦ Individual, swing away, removable legrest
- ♦ Grey molded plastic tray - 16" or 18" only
- ♦ Solid seat base with removable cushion
- ♦ 5", 6" or 8" heavy duty casters
- ♦ 6" or 9" brake extensions
- ♦ Side bolsters - 2"
- ♦ Pommel - velcro
- ♦ I.V. pole adjustable or fixed
- ♦ Oxygen tank holder

\*see full list of accessories on the MLTR 500 & 600 order form

## MODEL SPECIFICATIONS

- ♦ 5 year warranty on frame
- ♦ 2 year warranty on additional components



♦ TILT AND RECLINE



♦ FOLD DOWN BACK



♦ ADJUSTABLE LATERALS

## Dealer Information

6540 Tomken Road  
Mississauga ON L5T 2E9  
T 905.564.2250 F 905.564.5424  
1.800.362.0873

[www.mapleleafwheelchair.ca](http://www.mapleleafwheelchair.ca)

