

LECKEY[®]

Distributed by  Mobility

Advance Bath Chair



Secure Bathing Chair for Children, Teens & Young Adults.

Bathing is an important function in its own right, with the needs of the caregiver being as important as the needs of the child. The **Advance Bath Chair** provides users with appropriate levels of support to ensure safety and security during bathing.

Postural Support

Pelvic Stability

The **Advance Bath Chair** provides pelvic support with a soft, adjustable hip belt.



Promote Trunk Posture

Where additional trunk support is needed, a chest belt with integrated lateral supports can be attached.



Leg & Foot Positioning

Where legs have a tendency to scissor, the padded abduction belt can be attached, ensuring the legs stay in a symmetrical position.



Head Control

Head support can be easily added to the adjustable backrest to provide additional security during bathing.



Personalized Fit



The product can be set up into various positions depending on the individual's postural requirements, and the height can be adjusted to the appropriate height for the parent or caregiver to reduce the risk of back strain or injury.

Specifications

	Size 1	Size 2	Size 3	Size 4
Age	1-5	4-9	8-14	12-18
User Height	30"– 41"	37"– 49"	45"– 57"	53"– 65"
Seat Height	2" – 17"	2" – 17"	2" – 17"	2" – 17"
Weight Capacity	159 lbs.	159 lbs.	159 lbs.	159 lbs.
Total Width	17"	17"	17"	17"
Available Colors	Aquamarine, Royal Blue, Pink			

For complete specifications and available options, visit www.kimobility.com