

UNI-FIT EXTENDER

Increases SuperPole floor to ceiling range from 93"-99" to 100"-120".

Weight Capacity	300 lbs
Code	SP-UF

*Not to be used with SuperTrapeze, SuperPole HD, or SuperPole Ultra.




SUPERPOLE ULTRA

A longer version of the SuperPole with a floor to ceiling range from 121"-140".

Weight Capacity	300 lbs
Code	SP-ULTRA

*Not to be used with SuperBar, SuperTrapeze, Uni-Fit Extender, and Angled Ceiling Plate.



CEILING PLATE EXTENDER

Allows standard ceiling plate 20" length to span 30" with holes for ceiling attachment to 24" ceiling beam centers.

Code	SP-CPE
------	--------



ANGLED CEILING PLATE

Bolts to angled ceilings between 0-45 degrees. Top plate diameter is 18".

Code	SP-ACP (300lbs)
(Weight Capacity)	SuperPole with Angled Ceiling Plate
	SP-ACP-AK
	Angled Ceiling Plate Add-On Kit

*Not to be used with SuperTrapeze, SuperPole HD, or Ceiling Plate Extender.

