



# BRODA<sup>®</sup>

*Comfort is Essential*

**ELITE**  
TILT RECLINER

## INDUSTRY LEADING POSITIONING VERSATILITY AND FUNCTION

Enhanced positioning and innovative features provide exceptional comfort. Superior durability accommodates use in multiple care settings.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE (TRENDLENBURG CAPABLE)
- 3 HEIGHT ADJUSTABLE ARMRESTS
- 4 SWING AWAY REMOVABLE ARMS
- 5 REMOVABLE WINGS FOR UPPER LATERAL SUPPORT
- 6 SPLIT FLIP DOWN FOOTREST
- 7 MULTI-DIRECTIONAL LOCKING REAR CASTERS
- 8 INDEPENDENT ELEVATING, ARTICULATING, LENGTH ADJUSTABLE LEGREST

### PATIENT BENEFITS



Posterior tilt effectively opens the diaphragm improving digestion, oxygenation, and organ function.



Tilt in Space feature facilitates positioning for comfort, pain management and pain relief.



Lay flat feature facilitates supine transfers, reducing the risk associated with falls and injuries.



Provides ease of maneuverability to accommodate various activities (Middle Mag Wheel Option Shown).

### CAREGIVER AND FACILITY BENEFITS

- Facilitates positioning for comfort, pain management and improved respiratory function
- Optimum safety and functionality withstands severe agitation and involuntary movements
- Tilt and recline provide ease of use frequent re-positioning, helping improve patient outcomes
- Proper seating alignment decreases patient fatigue and postural deviations reducing the use of restraints
- Infinitely adjustable positioning cradles and soothes patient to reduce agitation and provide contentment

## COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue

## SPECIFICATIONS

|                      |                                                                                                                                |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Frame                | 16 Gauge tubular steel frame                                                                                                   |
| Seat & Backrest      | BRODA Comfort Tension Seating®                                                                                                 |
| Degrees of Tilt      | 40°                                                                                                                            |
| Degrees of Recline   | 90°(Trendelenburg capable)                                                                                                     |
| Mag Wheel Option     | 20" or 22" Front or rear mount (50cm, 55cm)<br>14" rear or mid mount (35cm)                                                    |
| Standard Seat Width  | 16", 18", 20" (40cm, 45cm, 50cm)                                                                                               |
| Custom Seat Width    | 22", 24", 26", 28" (55cm, 60cm, 65cm, 70cm)                                                                                    |
| Standard Seat Depth  | 17", 18.5" (43cm, 46cm)                                                                                                        |
| Custom Seat Depth    | 15.5", 20" (39cm, 50cm)                                                                                                        |
| Standard Seat Height | 20.5" (51cm)                                                                                                                   |
| Custom Seat Height   | 19", 22", 23.5" (48cm, 55cm, 59cm)                                                                                             |
| Front Casters        | 5" Front total lock (12.5cm)                                                                                                   |
| Rear Casters         | 6" Multi-directional locking (15cm)                                                                                            |
| Weight Capacity      | 350lbs. / 159kg                                                                                                                |
| Standard Padding     | Headrest cushion, shoulder bolsters, lower side pads,<br>upper side pads, seat and back pad with cover,<br>calf pad, sole pads |
| ADP Device Code*     | WMA5B-020, WMA5B-010                                                                                                           |

\* Specific configurations apply. Consult your rep for details

## COMMON ACCESSORIES

- Mag wheel package
- Aluminum front casters
- Twin wheel front casters
- Huntington's padding package (HSP)
- Padded armrests
- Full wrap vinyl seat and back covers
- BRODA app (additional positioning package)
- IV Pole
- O2 Holder
- ABS Legrest
- Solid or clear tray
- Half tray
- Padded tray cover
- Thigh belt



CTS Strapping flexes for individualized fit & superior comfort.

## UPHOLSTERY OPTIONS



RENAISSANCE  
BURGUNDY



RENAISSANCE  
CHARCOAL



RENAISSANCE  
METEOR



SILVERTEX  
SAPPHIRE BLUE

## FRAME OPTIONS



GREY  
HAMMER



BLACK



HEALTHCARE DIVISION

[info@brodaseating.com](mailto:info@brodaseating.com)

p. 1-844-552-7632

[Brodaseating.com](http://Brodaseating.com)

©Broda USA, INC.

8/16