

Bariatric Sentra EC Heavy-Duty Wheelchair Dual Cross Brace

HCPCS: K0007



- Weight Capacity: 450 lbs
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Double embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Standard dual axle provides easy transition of seat to hemi-level
- Urethane tires mounted on composite wheels provide durability and low maintenance
- Chrome hand rims standard
- 8" front casters, adjustable in three positions
- No maintenance, sealed, precision bearings
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD20ECDDAHD-SF	1/cs
STD20ECDDAHD-ELR	1/cs
STD20ECDFAH-D-SF	1/cs
STD20ECDFAH-D-ELR	1/cs
STD22ECDDA-SF	1/cs
STD22ECDDA-ELR	1/cs
STD22ECDF-A-SF	1/cs
STD22ECDF-A-ELR	1/cs
STD24ECDDA-SF	1/cs
STD24ECDDA-ELR	1/cs
STD24ECDF-A-SF	1/cs
STD24ECDF-A-ELR	1/cs

SPECIFICATIONS	Description
Width Open with 20" Seat	28"
Width Open with 22" Seat	30"
Width Open with 24" Seat	32"
Width Closed	12.5"
Height	36"
Seat-To-Floor Height	17.5" - 19.5"
Back Height (20" 22" and 24")	16" 18"
Seat Depth	18"
Length with Riggings	42"
Product Weight with 20" Seat	56 lbs.

SPECIFICATIONS	Description
Product Weight with 22" Seat	57 lbs.
Product Weight with 24" Seat	58 lbs.
Carton Shipping Weight with 20" Seat	62 lbs.
Carton Shipping Weight with 22" Seat	63 lbs.
Carton Shipping Weight with 24" Seat	64 lbs.
Weight Limit	450 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS801	Wheel Lock Extension, 1/ea
STDS827	Amputee Adapter, 1/pr
STDS806	Anti-Fold Bar, 1/ea
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS836	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS862	Universal Forward Stabilizer, 1/pr
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS801	Wheel Lock Extension, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS836	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS862	Universal Forward Stabilizer, 1/pr
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS1S2630	1/ea
STDS1S2632	1/ea
STDS1S2634	1/ea
STDS001071	
STDS1S1000HD	1/ea
STDS1S1012	1/ea