

The **Invacare® Top End® Crossfire™ T7A Wheelchair** features an adjustable, performance-minded, open frame design unlike any other rigid chair available today.

Features

Super light 7005 aluminum frame material
 Modular frame design allows for ample seat-to-floor height adjustments without re-squaring the casters
 16 degrees of back angle adjustments
 3.25" of center-of-gravity adjustments
 No charge mesh back upholstery offers more comfort in various climates
 Pressed bearing allows for a more controlled rotation and eliminates caster flutter

HCPCS Code: K0005



1. Adjustable angle fold down back
2. Invacare® Matrix® Seating available



Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Back Angle	Front Frame Angle	Camber	Tubing Diameter
CT7A	Seat Width + 6.75"	12"-20" W 14"-20" D	16"-21" Front 14"- 20" Rear .5" Increments	75°, 79°, 83°, 87°, 91°	70°, 75°, 80°, 85°	0°, 3°, 6°	1.12"
	Knee-to-Heel Range	Center-of-Gravity	Wheels	Casters	Total Weight**	Weight Capacity	Limited Warranty
	11.25"-20.25"	0.75"-3.75"	22", 24", 25", 26"	3", 4", 5"	17.8 lb.	250 lb.	5 Years Frame

**16"x16" with 24" lightweight performance spoke wheels, Primo Hi-Pressure tires, screw-on seat upholstery, mesh back upholstery with adjustable tension, long seat support and 4" urethane casters.

Please see the Owner's Manual or Invacare Configurator for complete option/specification combinations.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.